

# Acquisition and Analysis of Biosignals

## DTEK0042

Origins, Significance and Acquisition of Biosignals:  
Cardiovascular System

**Iman Azimi, Ph.D. (Tech.)**

Email: [iman.azimi@utu.fi](mailto:iman.azimi@utu.fi)



Turun yliopisto  
University of Turku


# Introduction

In the previous session, we learned:


- ❑ About the action potential, i.e., the origin of electric signal
- ❑ Different biosignals, e.g., EMG, EEG and EDA

In this session, we will learn:

- ❑ Biosignals collected to monitor the cardiovascular system
  - ECG, PCG, SCG, and PPG




# Electrocardiogram (ECG)



# Electrocardiogram (ECG)


- ❑ ECG is an electrical signal generated from the activities of heart
- ❑ ECG is the most commonly known & used biomedical signal
- ❑ The rhythm of the heart is estimated by counting the readily identifiable waves.
- ❑ ECG wave shape is altered by cardiovascular diseases and abnormalities (arrhythmias)



# The heart (1)


## ❑ A four-chambered pump:

1. **Right atrium (RA):** collect deoxygenated blood from the systemic circulation (body)
2. **Right ventricle (RV):** pump out to the lungs
3. **Left atrium (LA):** collect oxygenated blood from the lungs
4. **Left ventricle (LV):** pump out to the systemic circulation (body)



# The heart (2)


- ❑ In general there are two phases:
  1. **Systole**: contracting or pumping phase
  2. **Diastole**: resting or filling phase of a cardiac chamber
- ❑ Cardiac rhythm or Heart rate (HR) is controlled by the autonomic nervous system (ANS)
- ❑ **For the pumping**: There are two regions with specialized pacemaker cells:
  1. Sinoatrial (SA) node
  2. Atrioventricular (AV) node



# The electrical system of the heart

## ❑ In each cardiac cycle:

1. The SA node triggers its own train of action potentials
2. Electrical activity propagates through atria at comparatively low rates => slow-moving depolarization (contraction) of atria
3. The AV node fires (Propagation delay at the AV node => transfer blood from the atria to ventricle)
4. The bundle branches and Purkinje system propagate the stimulus
5. Electrical activity propagates through ventricles at a high rate => depolarization (contraction) of ventricle
6. A normally isoelectric segment is caused after the action potential of the ventricular muscle cells
7. Repolarization (relaxation) of the ventricles is occurred



# The ECG signal component

□ In a regular cardiac cycle, we measure:

**P-wave:** depolarization of atria

- (0.1 – 0.2 mV , 60 – 80 ms)

**PQ segment:** delay at the AV node

- (60 – 80 ms)

**QRS complex:** depolarization of ventricle


- (1 mV , 80 ms)

**ST segment:** a normally isoelectric segment

- (100 – 120 ms)

**T-wave:** repolarization of the ventricles

- (0.1 – 0.3 mV, 120 – 160 ms)







# Arrhythmia: disturbance in the cardiac cycle

- ❑ Irregular firing patterns from the SA node
- ❑ Abnormal and additional pacing activity from other parts of the heart
- ❑ If the SA node is depressed or inactive, other parts (e.g., AV node, purkinje fibers) may take over the role of the pacemaker
- ❑ Different arrhythmias result from variations in the site and frequency of impulse formation
  - Premature ventricular contractions => ventricular fibrillation
- ❑ Bundle-branch block causes a widened and jagged QRS
- ❑ ST segment might be elevated or depressed due to reduced blood supply to a part of the heart muscles (a block in the coronary arteries) or dead tissue

# Arrhythmia – example (1)



- ❑ Ventricular fibrillation in third and sixth beats



# Arrhythmia – example (2)

## ❑ Right bundle-branch block and hypertrophy


- QRS complex is wider than normal



# ECG signal acquisition


## □ Standard 12-channel ECG obtained

- Using 4 **limb leads**
  - Left arm
  - Right arm
  - Left leg
  - Right leg (ground)
- Using 6 **chest leads**



# ECG acquisition – limb (1)

- ❑ 6 channels are obtained using the 4 limb leads
- ❑ We assume a hypothetical equilateral triangle named as “**Einthoven’s triangle**”
- ❑ Center of the triangle is **Wilson’s central terminal** which is:
  - Formed by combining left arm, right arm, and left leg leads
  - Reference for chest leads
  - Schematically, location of the heart




Rangayyan, R. M. *Biomedical signal analysis*. 2nd Edition, Vol. 33. John Wiley & Sons, 2015.

# ECG acquisition – limb (2)

❑ **6 channels** are obtained using the 4 limb leads

- Lead I = LA – RA
- Lead II = LL – RA
- Lead III = LL – LA
- aVR
- aVL
- aVF


❑ The leads measure projections of the 3D cardiac electrical vector on to the axes of the leads.





[https://en.wikipedia.org/wiki/File:Limb\\_leads\\_of\\_EKG.png](https://en.wikipedia.org/wiki/File:Limb_leads_of_EKG.png)

# ECG acquisition – chest

- ☐ **6 channels** are obtained using the 6 chest leads
  - V1
  - V2
  - V3
  - V4
  - V5
  - V6
- ☐ V1 and V2: activity in the right-half of the heart
- ☐ V3 and V4: septal activity
- ☐ V5 and V6: left-ventricular activity



# 12-channel ECG - example






# Acquisition - notes


- ❑ In spite of being redundant, the 12-lead system serves as the basis of the standard clinical ECG.
- ❑ Clinical ECG for diagnosis
  - filtered to 0.05 – 100 Hz bandwidth
  - Sampling rate: 500 Hz
- ❑ ECG for heart-rate monitoring
  - filtered to 0.05 – 50 Hz bandwidth

# Phonocardiogram (PCG)




# Phonocardiogram (PCG)

- ❑ PCG is the sound signal generated from the activities of heart.
- ❑ The sound was traditionally collected by stethoscopes.
- ❑ Recording the PCG requires a transducer to convert the sound signal into an electronic signal => E.g., microphones and pressure transducers
- ❑ **Diagnosis:** Cardiovascular diseases cause changes or additional sounds and murmurs



# Heart sounds

- ❑ Heart sounds are the vibrations of the whole cardiovascular system triggered by pressure gradients.
- ❑ Heart sound components are best heard at certain locations on the chest
- ❑ 4 auscultatory areas:
  - Mitral area
  - Aortic area
  - Pulmonary area
  - Tricuspid area




# PCG signal acquisition

- ❑ Using a Piezoelectric contact sensors
  - bandwidth 0.05 – 1, 000 Hz
- ❑ PCG recording performed in a quiet room
- ❑ Patient is in supine position
- ❑ PCG transducer placed firmly at the desired position on the chest

# PCG signal components

- ❑ **Intensity, frequency content, and timing** of the heart sounds are important features
- ❑ A normal cardiac cycle contains two major sounds:
  - **S1**: Occurs at the onset of ventricular contraction
    - Caused by the closure of the heart valves, oscillation of the blood, etc.
  - **S2**: Closure of the semilunar valves (pulmonary and aortic valves)



Rangayyan, R. M. *Biomedical signal analysis*. 2nd Edition, Vol. 33. John Wiley & Sons, 2015.

# Other sounds

## ❑ Other sounds:


- S3: Sudden termination of the ventricular rapid-filling phase
- S4: atrial contractions displacing blood into the distended ventricles.
- Valvular clicks and snaps.

## ❑ Murmurs:

- Caused by cardiovascular diseases
- Are systolic and diastolic murmurs
- Are high-frequency, noise-like sounds
- Arise when the velocity of blood becomes high (flows through an irregularity)

# PCG and ECG


- ❑ ECG is electrical and PCG is mechanical (sound)
- ❑ Analysis of ECG is easier, especially if there are murmurs in the PCG.
- ❑ ECG can be used in the analysis of PCG
  - Use the QRS complex in the ECG to detect S1 in PCG



Rangayyan, R. M. *Biomedical signal analysis*. 2nd Edition, Vol. 33. John Wiley & Sons, 2015.




# Seismocardiogram (SCG)



# Seismocardiography (SCG)


- ❑ SCG is the signal which represents the local vibrations generated by the heart's contraction and the ejection of blood from the ventricles.
- ❑ SCG is characterized by several peaks and valleys reflecting specific events of the beating heart
- ❑ SCG can be used to detect heart rate and to diagnose diseases.



Lee, H., Lee, H., & Whang, M. (2018). An enhanced method to estimate heart rate from seismocardiography via ensemble averaging of body movements at six degrees of freedom. *Sensors*, 18(1), 238.

# SCG signal acquisition

- ❑ SCG can readily be detected by placing a low-noise accelerometer on the chest
  - Using wearables
- ❑ The precise location of the sensor on the chest impacts the measured signal
- ❑ A widely used position is on the sternum
- ❑ Patient is in supine position




Inan, Omer T., et al. "Novel wearable seismocardiography and machine learning algorithms can assess clinical status of heart failure patients." *Circulation: Heart Failure* 11.1 (2018): e004313.

# SCG signal acquisition

□ SCG is presented in three axes (i.e., X,Y, and Z)


- Most studies focus on **Z axis** (dorso-ventral)
- However, additional biological information could be derived also from the X axis (right-left) and Y axis (head-foot).



Lee, H., Lee, H., & Whang, M. (2018). An enhanced method to estimate heart rate from seismocardiography via ensemble averaging of body movements at six degrees of freedom. *Sensors*, 18(1), 238.


# SCG signal components

- ❑ **Amplitude, timing, or morphology** of the signal are important features
- ❑ A normal cardiac cycle contains peaks and valleys showing heart functions
  - E.g.: atrial systole (AS), mitral valve closure (MC), and aortic valve opening (AO)



# ECG, PCG, and SCG

- ❑ ECG is electrical, PCG is acoustic, and SCG is vibrational
- ❑ SCG is more sensitive to motion artifacts
- ❑ SCG shows the opening/closure of heart valves



Choudhary, T., Sharma, L. N., & Bhuyan, M. K. (2018). Heart sound extraction from sternal seismocardiographic signal. *IEEE Signal Processing Letters*, 25(4), 482-486.

# Relevant signals to SCG

## ❑ Gyrocardiogram<sup>1</sup>


- It shows heart motions using a sensor of angular motion, **gyroscope**, attached to the skin of the chest.

## ❑ Ballistocardiogram<sup>2</sup>

- It shows the movement of the body caused by cardiac ejection of blood using an **accelerometer** sensor attached under the body (on the bed).

1. M. Jafari Tadi, et al. "Gyrocardiography: A new non-invasive monitoring method for the assessment of cardiac mechanics and the estimation of hemodynamic variables." Scientific reports 7.1 (2017): 6823.
2. Omer T. Inan, et al. "Ballistocardiography and seismocardiography: A review of recent advances." IEEE journal of biomedical and health informatics 19.4 (2014): 1414-1427.


# Photoplethysmogram (PPG)






# Photoplethysmogram (PPG)

- ❑ PPG is the signal which indicates blood volumetric changes in the tissue
- ❑ The optical characteristics of red blood cells changes according to the amount of oxygen it carries via hemoglobin



# PPG signal acquisition

- ❑ **Emitter or light source:** to emit light to the skin
- ❑ **Detector or light sensor:** to absorb the reflected light
- ❑ The hemoglobin light absorption coefficient differs with different light wavelength:
  - Red
  - Infrared
  - Green



# Wrist-based PPG sensor


- ❑ One light source: **green**
- ❑ Green LED has much greater absorptivity for both oxyhemoglobin and deoxyhemoglobin
  - Better signal-to-noise ratio
- ❑ Used in wristbands and smart watches
- ❑ Heart rate monitoring



<https://buy.garmin.com/en-US/US/p/567813/pn/010-01755-10#overview>


# Finger-based PPG sensor

- ❑ Two light sources: **red** and **infrared**
- ❑ Two PPG signals are collected with different intensity due to the hemoglobin light absorption coefficient
  - $SpO_2$  or blood oxygen saturation can be obtained.



Amiri, D., et al. "Optimizing energy in wearable devices using fog computing." Fog Computing: Theory and Practice, (2019).

# Vital signs detection




Kasaeyan Naeini, E., et al. "A Real-time PPG Quality Assessment Approach for Healthcare Internet-of-Things." *Procedia Computer Science* 151 (2019): 551-558.

- ❑ Variations in the PPG signal are associated with the oscillation of heartbeat and respiration
  - Heart rate, heart rate variability, and respiration rate can be tracked
- ❑ Finger-based sensor:  $\text{SpO}_2$
- ❑ Wrist-based sensor: signal-to-noise ratio is too low to detect respiration rate


# ECG and PPG

- ❑ ECG is electrical and PPG is light reflection.
- ❑ Analysis of ECG and PPG may provide new information
- ❑ E.g., Pulse transit time (PTT)
  - Interval between the peak of the R-wave in ECG and the peak in PPG
  - Can be used to estimate systolic and diastolic blood pressure



Ma, H. T. (2014). A blood pressure monitoring method for stroke management. BioMed research international, 2014.

Wang, R., et al. (2014, October). Cuff-free blood pressure estimation using pulse transit time and heart rate. In 2014 12th international conference on signal processing (ICSP) (pp. 115-118). IEEE.



# Conclusion

In this session, we learned:

- ☐ Origin and acquisition of biosignals, collected to monitor the cardiovascular system

In the next session, we will learn:

- ☐ About noise and artifacts in biosignals
- ☐ Filtering techniques



# Thank You

## Questions?



Turun yliopisto  
University of Turku